

THADA Antibody



PACO44529

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

THADA

Uniprot

Q6YHU6

Synonyms:

Thyroid adenoma-associated protein (Gene inducing thyroid adenomas protein), THADA, GITA KIAA1767


Immunogen:

Recombinant Human Thyroid adenoma-associated protein (1684-1953AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO44529 at dilution of 1:100.