
Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Asymmetrically hydrolyzes Ap4A to yield AMP and ATP. Plays a major role in maintaining homeostasis.

Gene ID:

NUDT2

Uniprot

P50583

Synonyms:

Bis(5'-nucleosyl)-tetraphosphatase [asymmetrical] (EC 3.6.1.17) (Diadenosine 5',5'''-P1, P4-tetraphosphate asymmetrical hydrolase) (Ap4A hydrolase) (Ap4Aase) (Diadenosine tetraphosphatase) (Nucleoside diphosphate-linked moiety X motif 2) (Nudix motif 2), NUDT2, APAH1


Immunogen:

Recombinant Human Bis(5'-nucleosyl)-tetraphosphatase [asymmetrical] protein (1-147AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO44545 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human breast cancer using PACO44545 at dilution of 1:100.