

TP53AIP1 Antibody



PACO44546

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May play an important role in mediating p53/TP53-dependent apoptosis.

Gene ID:

TP53AIP1

Uniprot

Q9HCN2

Synonyms:

p53-regulated apoptosis-inducing protein 1 (p53AIP1), TP53AIP1


Immunogen:

Recombinant Human p53-regulated apoptosis-inducing protein 1 protein (1-108AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human liver cancer using PACO44546 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human thymus tissue using PACO44546 at dilution of 1:100.