

PACO44548

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Participates in chain elongation of fatty acid, . Has no 2,4-dienoyl-CoA reductase activity.

Gene ID:

PECR

Uniprot

Q9BY49

Synonyms:

Peroxisomal trans-2-enoyl-CoA reductase (TERP) (EC 1.3.1.38) (2,4-dienoyl-CoA reductase-related protein) (DCR-RP) (HPDHase) (Short chain dehydrogenase/reductase family 29C member 1) (pVI-ARL), PECR, SDR29C1


Immunogen:

Recombinant Human Peroxisomal trans-2-enoyl-CoA reductase protein (24-303AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO44548 at dilution of 1:100.