

TXNL4B Antibody



PACO44563

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Essential role in pre-mRNA splicing. Required in cell cycle progression for S/G(2) transition.

Gene ID:

TXNL4B

Uniprot

Q9NX01

Synonyms:

Thioredoxin-like protein 4B (Dim1-like protein), TXNL4B, DIM2 DLP


Immunogen:

Recombinant Human Thioredoxin-like protein 4B protein (1-149AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO44563 at dilution of 1:100.