

BLID Antibody



PACO44568

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Functions as a proapoptotic molecule through the caspase-dependent mitochondrial pathway of cell death.

Gene ID:

BLID

Uniprot

Q8IZY5

Synonyms:

BH3-like motif-containing cell death inducer (Breast cancer cell protein 2), BLID, BRCC2


Immunogen:

Recombinant Human BH3-like motif-containing cell death inducer protein (1-108AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO44568 at dilution of 1:100.