

GCM1 Antibody



PACO44579

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Transcription factor that is necessary for placental development. Binds to the trophoblast-specific element 2 (TSE2) of the aromatase gene enhancer.

Gene ID:

GCM1

Uniprot

Q9NP62

Synonyms:

Chorion-specific transcription factor GCMA (hGCMA) (GCM motif protein 1) (Glial cells missing homolog 1), GCM1, GCMA

Immunogen:

Recombinant Human Chorion-specific transcription factor GCMA protein (207-436AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot. All lanes: GCM1 antibody at 0.97 µg/ml. Lane 1: HeLa whole cell lysate. Lane 2: HepG2 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 49 kDa. Observed band size: 49 kDa.