

LMO4 Antibody



PACO44587

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Probable transcriptional factor.

Gene ID:

LMO4

Uniprot

P61968

Synonyms:

LIM domain transcription factor LMO4 (Breast tumor autoantigen) (LIM domain only protein 4) (LMO-4), LMO4


Immunogen:

Recombinant Human LIM domain transcription factor LMO4 protein (1-165AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO44587 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO44587 at dilution of 1:100.