

PDCD2L Antibody



PACO44592

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:

Over-expression suppresses AP1, CREB, NFAT, and NF-kB transcriptional activation, and delays cell cycle progression at S phase.

Gene ID:

PDCD2L

Uniprot

Q9BRP1

Synonyms:

Programmed cell death protein 2-like, PDCD2L


Immunogen:

Recombinant Human Programmed cell death protein 2-like protein (1-358AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot. All lanes: PDCD2L antibody at 1.69 μ g/ml. Lane 1: Hela whole cell lysate. Lane 2: MCF-7 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 39 kDa. Observed band size: 39 kDa.



Immunohistochemistry of paraffin-embedded human liver cancer using PACO44592 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO44592 at dilution of 1:100.