

PACO44604

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

E3 ubiquitin-protein ligase mediating 'Lys-48'-and 'Lys-63'-linked polyubiquitination and subsequent proteasomal degradation of modified proteins. May play a role in myelopoiesis.

Gene ID:

ARIH2

Uniprot

O95376

Synonyms:

E3 ubiquitin-protein ligase ARIH2 (ARI-2) (Protein ariadne-2 homolog) (EC 2.3.2.27) (RING-type E3 ubiquitin transferase ARIH2) (Triad1 protein), ARIH2, ARI2 TRIAD1


Immunogen:

Recombinant Human E3 ubiquitin-protein ligase ARIH2 protein (194-493AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human lung tissue using PACO44604 at dilution of 1:100.