

ATOH7 Antibody



PACO44607

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Transcription factor involved in the differentiation of retinal ganglion cells.

Gene ID:

ATOH7

Uniprot

Q8N100

Synonyms:

Protein atonal homolog 7 (Class A basic helix-loop-helix protein 13) (bHLHa13) (Helix-loop-helix protein hATH-5) (hATH5), ATOH7, ATH5 BHLHA13


Immunogen:

Recombinant Human Protein atonal homolog 7 protein (47-152AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO44607 at dilution of 1:100.