

# ATOH7 Antibody



PACO44608

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Transcription factor involved in the differentiation of retinal ganglion cells.

**Gene ID:**

ATOH7

**Uniprot**

Q8N100

**Synonyms:**

Protein atonal homolog 7 (Class A basic helix-loop-helix protein 13) (bHLHa13) (Helix-loop-helix protein hATH-5) (hATH5), ATOH7, ATH5 BHLHA13

**Immunogen:**


Recombinant Human Protein atonal homolog 7 protein (47-152AA).

**Storage:**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human colon cancer using PACO44608 at dilution of 1:100.