

BNIP1 Antibody



PACO44611

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

SNARE that may be involved in targeting and fusion of Golgi-derived retrograde transport vesicles with the ER. Required for maintenance of ER network. Implicated in the suppression of cell death. May be involved in mitochondrial autophagy.

Gene ID:

BNIP1

Uniprot

Q12981

Synonyms:

Vesicle transport protein SEC20 (BCL2/adenovirus E1B 19 kDa protein-interacting protein 1) (Transformation-related gene 8 protein) (TRG-8), BNIP1, NIP1 SEC20L


Immunogen:

Recombinant Human Vesicle transport protein SEC20 protein (1-200AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human liver tissue using PACO44611 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO44611 at dilution of 1:100.