

DNAJB9 Antibody



PACO44622

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Involved in endoplasmic reticulum-associated degradation (ERAD) of misfolded proteins. Acts as a co-chaperone with an Hsp70 protein.

Gene ID:

DNAJB9

Uniprot

Q9UBS3

Synonyms:

DnaJ homolog subfamily B member 9 (Endoplasmic reticulum DNA J domain-containing protein 4) (ER-resident protein ERdj4) (ERdj4) (Microvascular endothelial differentiation gene 1 protein) (Mdg-1), DNAJB9, MDG1


Immunogen:

Recombinant Human DnaJ homolog subfamily B member 9 protein (24-223AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO44622 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO44622 at dilution of 1:100.