

PEX3 Antibody



PACO44671

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:

Involved in peroxisome biosynthesis and integrity. Assembles membrane vesicles before the matrix proteins are translocated. As a docking factor for PEX19, is necessary for the import of peroxisomal membrane proteins in the peroxisomes.

Gene ID:

PEX3

Uniprot

P56589

Synonyms:

Peroxisomal biogenesis factor 3 (Peroxin-3) (Peroxisomal assembly protein PEX3), PEX3


Immunogen:

Recombinant Human Peroxisomal biogenesis factor 3 protein (144-373AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot. All lanes: PEX3 antibody at 1.24 μ g/ml + Mouse kidney tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 42 kDa. Observed band size: 42 kDa.



Immunohistochemistry of paraffin-embedded human liver tissue using PACO44671 at dilution of 1:100.