

SUPT4H1 Antibody



PACO44696

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

SUPT4H1

Uniprot

P63272

Synonyms:

Transcription elongation factor SPT4 (hSPT4) (DRB sensitivity-inducing factor 14 kDa subunit) (DSIF p14) (DRB sensitivity-inducing factor small subunit) (DSIF small subunit), SUPT4H1, SPT4H SUPT4H


Immunogen:

Recombinant Human Transcription elongation factor SPT4 protein (1-117AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO44696 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO44696 at dilution of 1:100.