

TRIM29 Antibody



PACO44705

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

It is able to complement the radiosensitivity defect of an ataxia telangiectasia (AT) fibroblast cell line.

Gene ID:

TRIM29

Uniprot

Q14134

Synonyms:

Tripartite motif-containing protein 29 (Ataxia telangiectasia group D-associated protein), TRIM29, ATDC


Immunogen:

Recombinant Human Tripartite motif-containing protein 29 protein (309-588AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human skin tissue using PACO44705 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO44705 at dilution of 1:100.