

PACO44711

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

UDPGT is of major importance in the conjugation and subsequent elimination of potentially toxic xenobiotics and endogenous compounds.

**Gene ID:**

UGT2B10

**Uniprot**

P36537

**Synonyms:**

UDP-glucuronosyltransferase 2B10 (UDPGT 2B10) (EC 2.4.1.17), UGT2B10

**Immunogen:**


Recombinant Human UDP-glucuronosyltransferase 2B10 protein (20-280AA).

**Storage:**


PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human cervical cancer using PACO44711 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human ovarian cancer using PACO44711 at dilution of 1:100.