

ZBTB25 Antibody



PACO44725

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZBTB25

Uniprot

P24278

Synonyms:

Zinc finger and BTB domain-containing protein 25 (Zinc finger protein 46) (Zinc finger protein KUP), ZBTB25, C14orf51 KUP ZNF46


Immunogen:

Recombinant Human Zinc finger and BTB domain-containing protein 25 protein (266-435AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot. All lanes: ZBTB25 antibody at 0.43 μ g/ml. Lane 1: MCF-7 whole cell lysate. Lane 2: Hela whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 49 kDa. Observed band size: 49 kDa.



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO44725 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO44725 at dilution of 1:100.