

ZNF169 Antibody



PACO44732

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF169

Uniprot

Q14929

Synonyms:

Zinc finger protein 169, ZNF169

Immunogen:

Recombinant Human Zinc finger protein 169 protein (1-85AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot. All lanes: ZNF169 antibody at 5.32 μ g/ml + HeLa whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 68 kDa. Observed band size: 68 kDa.



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO44732 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO44732 at dilution of 1:100.