

ZNF324 Antibody



PACO44742

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF324

Uniprot

O75467

Synonyms:

Zinc finger protein 324A (Zinc finger protein ZF5128), ZNF324, ZNF324A


Immunogen:

Recombinant Human Zinc finger protein 324A protein (1-270AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot. All lanes: ZNF324 antibody at 1.57 μ g/ml. Lane 1: A549 whole cell lysate. Lane 2: 293T whole cell lysate. Lane 3: MCF-7 whole cell lysate. Lane 4: Jurkat whole cell lysate. Secondary: Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 61 kDa. Observed band size: 61 kDa.



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO44742 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO44742 at dilution of 1:100.