

ZNF350 Antibody



PACO44743

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Transcriptional repressor. Binds to a specific sequence, 5'-GGGxxxCAGxxxTTT-3', within GADD45 intron 3.

Gene ID:

ZNF350

Uniprot

Q9GZX5

Synonyms:

Zinc finger protein 350 (KRAB zinc finger protein ZFQR) (Zinc finger and BRCA1-interacting protein with a KRAB domain 1) (Zinc finger protein ZBRK1), ZNF350, ZBRK1


Immunogen:

Recombinant Human Zinc finger protein 350 protein (233-532AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO44743 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO44743 at dilution of 1:100.