

PACO44752

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

CHMP4B

Uniprot

Q9H444

Synonyms:

Charged multivesicular body protein 4b (Chromatin-modifying protein 4b) (CHMP4b) (SNF7 homolog associated with Alix 1) (SNF7-2) (hSnf7-2) (Vacuolar protein sorting-associated protein 32-2) (Vps32-2) (hVps32-2), CHMP4B, C20orf178 SHAX1


Immunogen:

Recombinant Human Charged multivesicular body protein 4b protein (1-224AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO44752 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human heart tissue using PACO44752 at dilution of 1:100.