
Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IP

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200, IP:1:200-1:2000

Protein Background:

Catalyzes the phosphorylation of the ketose sugar fructose to fructose-1-phosphate.

Gene ID:

KHK

Uniprot

P50053

Synonyms:

Ketohexokinase (EC 2.7.1.3) (Hepatic fructokinase), KHK


Immunogen:

Recombinant Human Ketohexokinase protein (1-298AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot. All lanes: KHKantibody at 3.6 μ g/ml. Lane 1: Mouse liver tissue. Lane 2: Mouse kidney tissue. Lane 3: Rat kidney tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 33 kDa. Observed band size: 33 kDa.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO44773 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human liver tissue using PACO44773 at dilution of 1:100.