
Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

RNA-binding protein that plays an essential role in spermatogenesis. May act by binding to the 3'-UTR of mRNAs and regulating their translation.

Gene ID:

DAZ2

Uniprot

Q13117

Synonyms:

Deleted in azoospermia protein 2, DAZ2


Immunogen:

Recombinant Human Deleted in azoospermia protein 2 protein (1-200AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO44826 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human skin tissue using PACO44826 at dilution of 1:100.