

ICAM4 Antibody



PACO44831

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

ICAM proteins are ligands for the leukocyte adhesion protein LFA-1 (integrin alpha-L/beta-2). ICAM4 is also a ligand for alpha-4/beta-1 and alpha-V integrins.

Gene ID:

ICAM4

Uniprot

Q14773

Synonyms:

Intercellular adhesion molecule 4 (ICAM-4) (Landsteiner-Wiener blood group glycoprotein) (LW blood group protein) (CD antigen CD242), ICAM4, LW


Immunogen:

Recombinant Human Intercellular adhesion molecule 4 protein (20-140AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human liver cancer using PACO44831 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO44831 at dilution of 1:100.