

PLA2G4C Antibody



PACO44845

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:

Has a preference for arachidonic acid, at the sn-2 position of phosphatidylcholine as compared with palmitic acid, .

Gene ID:

PLA2G4C

Uniprot

Q9UP65

Synonyms:

Cytosolic phospholipase A2 gamma (cPLA2-gamma) (EC 3.1.1.4) (Phospholipase A2 group IVC), PLA2G4C


Immunogen:

Recombinant Human Cytosolic phospholipase A2 γ protein (1-189AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot. All lanes: COX10 antibody at 5.84µg/ml + Mouse heart tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 61, 60, 63 kDa. Observed band size: 61 kDa.



Immunohistochemistry of paraffin-embedded human heart tissue using PACO44845 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO44845 at dilution of 1:100.