

COX10 Antibody



PACO44868

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Converts protoheme IX and farnesyl diphosphate to heme O.

Gene ID:

COX10

Uniprot

Q12887

Synonyms:

Protoheme IX farnesyltransferase, mitochondrial (EC 2.5.1. -) (Heme O synthase), COX10


Immunogen:

Recombinant Human Protoheme IX farnesyltransferase, mitochondrial protein (1-160AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO44868 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO44868 at dilution of 1:100.