

ANKRD1 Antibody



PACO44918

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:

May play an important role in endothelial cell activation. May act as a nuclear transcription factor that negatively regulates the expression of cardiac genes. Induction seems to be correlated with apoptotic cell death in hepatoma cells.

Gene ID:

ANKRD1

Uniprot

Q15327

Synonyms:

Ankyrin repeat domain-containing protein 1 (Cardiac ankyrin repeat protein) (Cytokine-inducible gene C-193 protein) (Cytokine-inducible nuclear protein), ANKRD1, C193 CARP HA1A2


Immunogen:

Recombinant Human Ankyrin repeat domain-containing protein 1 protein (1-130AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot. All lanes: ANKRD1 antibody at 1.91 μ g/ml. Lane 1: HeLa whole cell lysate. Lane 2: 293T whole cell lysate. Lane 3: Jurkat whole cell lysate. Lane 4: HepG2 whole cell lysate. Secondary: Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 36 kDa. Observed band size: 36 kDa.



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO44918 at dilution of 1:100.