

FMO5 Antibody



PACO44939

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:

In contrast with other forms of FMO it does not seem to be a drug-metabolizing enzyme.

Gene ID:

FMO5

Uniprot

P49326

Synonyms:

Dimethylaniline monooxygenase [N-oxide-forming] 5 (EC 1.14.13.8) (Dimethylaniline oxidase 5) (Hepatic flavin-containing monooxygenase 5) (FMO 5), FMO5


Immunogen:

Recombinant Human Dimethylaniline monooxygenase [N-oxide-forming] 5 protein (1-285AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot. All lanes: FMO5 antibody at 6.04 μ g/ml + Mouse liver tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 61, 33, 53 kDa. Observed band size: 61 kDa.



Immunohistochemistry of paraffin-embedded human liver tissue using PACO44939 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human lung tissue using PACO44939 at dilution of 1:100.