

# DCTN3 Antibody



PACO45070

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Together with dynein may be involved in spindle assembly and cytokinesis.

**Gene ID:**

DCTN3

**Uniprot**

O75935

**Synonyms:**

Dynactin subunit 3 (Dynactin complex subunit 22 kDa subunit) (p22), DCTN3, DCTN22

**Immunogen:**


Recombinant Human Dynactin subunit 3 protein (1-176AA).

**Storage:**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO45070 at dilution of 1:100.