
Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

As part of the NSL complex it is involved in acetylation of nucleosomal histone H4 on several lysine residues and therefore may be involved in the regulation of transcription.

Gene ID:

KANSL3

Uniprot

Q9P2N6

Synonyms:

KAT8 regulatory NSL complex subunit 3 (NSL complex protein NSL3) (Non-specific lethal 3 homolog) (Serum inhibited-related protein) (Testis development protein PRTD), KANSL3, KIAA1310 NSL3 PRTD SI1


Immunogen:

Recombinant Human KAT8 regulatory NSL complex subunit 3 protein (1-200AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO45091 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO45091 at dilution of 1:100.