

NEDD1 Antibody



PACO45099

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IP

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200, IP:1:200-1:2000

Protein Background:

Required for mitosis progression. Promotes the nucleation of microtubules from the spindle.

Gene ID:

NEDD1

Uniprot

Q8NHV4

Synonyms:

Protein NEDD1 (Neural precursor cell expressed developmentally down-regulated protein 1) (NEDD-1), NEDD1


Immunogen:

Recombinant Human Protein NEDD1 protein (411-660AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot. All lanes: NEDD1 antibody at 3.25µg/ml. Lane 1: Mouse small intestine tissue. Lane 2: Jurkat whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 72, 63, 73 kDa. Observed band size: 72 kDa.



Immunohistochemistry of paraffin-embedded human thyroid tissue using PACO45099 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human glioma using PACO45099 at dilution of 1:100.