
Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

cytosol, fructokinase activity, glucokinase activity, hexokinase activity, mannokinase activity, canonical glycolysis, cellular glucose homeostasis, glucose transport, glycolytic process.

Gene ID:

HK3

Uniprot

P52790

Synonyms:

Hexokinase-3 (EC 2.7.1.1) (Hexokinase type III) (HK III), HK3


Immunogen:

Recombinant Human Hexokinase-3 protein (1-270AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human lung cancer using PACO45141 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human spleen tissue using PACO45141 at dilution of 1:100.