
Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:

Substrate-recognition component of the SCF (SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex. Promotes the proteasome-dependent degradation of key sarcomeric proteins, such as alpha-actinin (ACTN2) and filamin-C (FLNC), essential for maintenance of normal contractile function.

Gene ID:

FBXO22

Uniprot

Q8NEZ5

Synonyms:

F-box only protein 22 (F-box protein FBX22p44), FBXO22, FBX22


Immunogen:

Recombinant Human F-box only protein 22 protein (104-403AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot. All lanes: FBXO22 antibody at 1.98 μ g/ml. Lane 1: Mouse heart tissue. Lane 2: Hela whole cell lysate. Lane 3: MCF-7 whole cell lysate. Lane 4: HepG2 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 45, 5, 31 kDa. Observed band size: 45 kDa.



Immunohistochemistry of paraffin-embedded human liver cancer using PACO45181 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO45181 at dilution of 1:100.