

# RNF13 Antibody



PACO45198

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

E3 ubiquitin-protein ligase that may play a role in controlling cell proliferation.

**Gene ID:**

RNF13

**Uniprot**

O43567

**Synonyms:**

E3 ubiquitin-protein ligase RNF13 (EC 2.3.2.27) (RING finger protein 13) (RING-type E3 ubiquitin transferase RNF13), RNF13, RZF

**Immunogen:**


Recombinant Human E3 ubiquitin-protein ligase RNF13 protein (204-381AA).

**Storage:**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO45198 at dilution of 1:100.