

RHOBTB1 Antibody



PACO45265

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

cytosol, plasma membrane, regulation of small GTPase mediated signal transduction.

Gene ID:

RHOBTB1

Uniprot

O94844

Synonyms:

Rho-related BTB domain-containing protein 1, RHOBTB1, KIAA0740


Immunogen:

Recombinant Human Rho-related BTB domain-containing protein 1 protein (240-400AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human thyroid tissue using PACO45265 at dilution of 1:100.