

SAV1 Antibody



PACO45270

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000,
IHC:1:20-1:500, IF:1:50-1:200

Protein Background:**Gene ID:**

SAV1

Uniprot

Q9H4B6

Synonyms:

Protein salvador homolog 1 (45 kDa WW domain protein) (hWW45), SAV1, WW45


Immunogen:

Recombinant Human Protein salvador homolog 1 protein (1-210AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO45270 at dilution of 1:100.



Immunofluorescent analysis of A549 cells using PACO45270 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO45270 at dilution of 1:100.