

TRIP12 Antibody




PACO45278


Product Information

Size:	Protein Background:
50ul	
Reactivity:	Gene ID:
Human	TRIP12
Source:	Uniprot
Rabbit	Q14669
Isotype:	Synonyms:
IgG	E3 ubiquitin-protein ligase TRIP12 (EC 2.3.2.26) (E3 ubiquitin-protein ligase for Arf) (ULF) (HECT-type E3 ubiquitin transferase TRIP12) (Thyroid receptor-interacting protein 12) (TR-interacting protein 12) (TRIP-12), TRIP12, KIAA0045 ULF
Applications:	Immunogen:
ELISA, IHC	Recombinant Human E3 ubiquitin-protein ligase TRIP12 protein (1-240AA).
Recommended dilutions:	Storage:
ELISA:1:2000-1:10000, IHC:1:20-1:200	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO45278 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human melanoma using PACO45278 at dilution of 1:100.