

DUSP16 Antibody



PACO45288

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Dual specificity protein phosphatase involved in the inactivation of MAP kinases. Dephosphorylates MAPK10 bound to ARRB2.

Gene ID:

DUSP16

Uniprot

Q9BY84

Synonyms:

Dual specificity protein phosphatase 16 (EC 3.1.3.16) (EC 3.1.3.48) (Mitogen-activated protein kinase phosphatase 7) (MAP kinase phosphatase 7) (MKP-7), DUSP16, KIAA1700 MKP7


Immunogen:

Recombinant Human Dual specificity protein phosphatase 16 protein (416-665AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human lung tissue using PACO45288 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO45288 at dilution of 1:100.