

TRPM4 Antibody



PACO45305

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

TRPM4

Uniprot

Q8TD43

Synonyms:

Transient receptor potential cation channel subfamily M member 4 (hTRPM4) (Calcium-activated non-selective cation channel 1) (Long transient receptor potential channel 4) (LTrpC-4) (LTrpC4) (Melastatin-4), TRPM4, LTRPC4


Immunogen:

Recombinant Human Transient receptor potential cation channel subfamily M member 4 protein (1-240AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO45305 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO45305 at dilution of 1:100.