

ZIC3 Antibody



PACO45306

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Acts as transcriptional activator. Required in the earliest stages in both axial midline development and left-right (LR) asymmetry specification. Binds to the minimal GLI-consensus sequence 5'-GGGTGGTC-3'.

Gene ID:

ZIC3

Uniprot

O60481

Synonyms:

Zinc finger protein ZIC 3 (Zinc finger protein 203) (Zinc finger protein of the cerebellum 3), ZIC3, ZNF203


Immunogen:

Recombinant Human Zinc finger protein ZIC 3 protein (1-240AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO45306 at dilution of 1:100.