

# CORO1C Antibody



PACO45316

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

May be involved in cytokinesis, motility, and signal transduction.

**Gene ID:**

CORO1C

**Uniprot**

Q9ULV4

**Synonyms:**

Coronin-1C (Coronin-3) (hCRNN4), CORO1C, CRN2 CRNN4

**Immunogen:**


Recombinant Human Coronin-1C protein (375-474AA).

**Storage:**


PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human testis tissue using PACO45316 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human cervical cancer using PACO45316 at dilution of 1:100.