

CD8B Antibody



PACO45341

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Identifies cytotoxic/suppressor T-cells that interact with MHC class I bearing targets. CD8 is thought to play a role in the process of T-cell mediated killing.

Gene ID:

CD8B

Uniprot

P10966

Synonyms:

T-cell surface glycoprotein CD8 beta chain (CD antigen CD8b), CD8B, CD8B1


Immunogen:

Recombinant Human T-cell surface glycoprotein CD8 β chain protein (22-170AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO45341 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO45341 at dilution of 1:100.