

CLEC1B Antibody



PACO45344

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

integral component of plasma membrane, plasma membrane, transmembrane signaling receptor activity, cell surface receptor signaling pathway, defense response, platelet activation.

Gene ID:

CLEC1B

Uniprot

Q9P126

Synonyms:

C-type lectin domain family 1 member B (C-type lectin-like receptor 2) (CLEC-2), CLEC1B, CLEC2


Immunogen:

Recombinant Human C-type lectin domain family 1 member B protein (55-229AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human liver tissue using PACO45344 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human liver cancer using PACO45344 at dilution of 1:100.