

GABRA4 Antibody



PACO45345

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

GABA, the major inhibitory neurotransmitter in the vertebrate brain, mediates neuronal inhibition by binding to the GABA/benzodiazepine receptor and opening an integral chloride channel.

Gene ID:

GABRA4

Uniprot

P48169

Synonyms:

Gamma-aminobutyric acid, receptor subunit alpha-4 (GABA(A) receptor subunit alpha-4), GABRA4

Immunogen:

Recombinant Human Gamma-aminobutyric acid, receptor subunit alpha-4 protein (36-258AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO45345 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human brain tissue using PACO45345 at dilution of 1:100.