
Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

This is a receptor for GIP. The activity of this receptor is mediated by G proteins which activate adenylyl cyclase.

Gene ID:

GIPR

Uniprot

P48546

Synonyms:

Gastric inhibitory polypeptide receptor (GIP-R) (Glucose-dependent insulinotropic polypeptide receptor), GIPR


Immunogen:

Recombinant Human Gastric inhibitory polypeptide receptor protein (22-138AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO45388 at dilution of 1:100.



Immunofluorescent analysis of 293 cells using PACO45388 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO45388 at dilution of 1:100.