

FZD3 Antibody



PACO45439

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

FZD3

Uniprot

Q9NPG1

Synonyms:

Frizzled-3 (Fz-3) (hFz3), FZD3


Immunogen:

Recombinant Human Frizzled-3 protein (23-205AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO45439 at dilution of 1:100.