

# IL17RB Antibody



PACO45448

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Receptor for the proinflammatory cytokines IL17B and IL17E. May play a role in controlling the growth and/or differentiation of hematopoietic cells.

**Gene ID:**

IL17RB

**Uniprot**

Q9NRM6

**Synonyms:**

Interleukin-17 receptor B (IL-17 receptor B) (IL-17RB) (Cytokine receptor-like 4) (IL-17 receptor homolog 1) (IL-17Rh1) (IL17Rh1) (Interleukin-17B receptor) (IL-17B receptor), IL17RB, CRL4 EVI27 IL17BR

**Immunogen:**


Recombinant Human Interleukin-17 receptor B protein (18-292AA).

**Storage:**


PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human liver tissue using PACO45448 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO45448 at dilution of 1:100.