

NRXN1 Antibody



PACO45462

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

NRXN1

Uniprot

Q9ULB1

Synonyms:

Neurexin-1 (Neurexin I-alpha) (Neurexin-1-alpha), NRXN1, KIAA0578


Immunogen:

Recombinant Human Neurexin-1 protein (31-310AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO45462 at dilution of 1:100.